**Bovine Spongiform Encephalopathy-Scrapie  
Antigen Test Kit, EIA**

**Kit Encéphalopathie Spongiforme Bovine–Tremblante  
(Scrapie) Antigène, EIA**

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English Version

## Overview

The IDEXX HerdChek\* Bovine Spongiform Encephalopathy (BSE)-Scrapie Antigen Test Kit is an antigen-capture enzyme immunoassay (EIA) for detection of the abnormal conformer of the prion protein (PrP<sup>Sc</sup>) in postmortem brain (obex preferred) tissues from bovines, small ruminants (sheep and goats) and cervids affected by BSE, scrapie, or chronic wasting disease (CWD). For small-ruminant applications, the test will also detect PrP<sup>Sc</sup> in lymph node or spleen tissue; for cervids the test will also detect PrP<sup>Sc</sup> in retropharyngeal lymph nodes (RPLN). It is designed to rapidly identify samples containing disease-associated PrP<sup>Sc</sup> with minimal sample handling, and can be automated for high throughput applications. This kit is for veterinary use only.

For European Member States: This is an EU-approved rapid test for in vitro detection of BSE and scrapie-related PrP<sup>Sc</sup> in cattle and small ruminants. This test can also be used in the EU for rapid testing of cervids for transmissible spongiform encephalopathy (TSE). The validation of this kit was carried out using tonsil, retropharyngeal lymphoid node and brainstem samples from elk (*Cervus canadensis*).

The producer of the rapid test must have a quality assurance system in place agreed to by the European Union Reference Laboratory (EURL), which ensures that the test performance does not change. The producer must provide the test protocol to the EURL. Sampling tools and modifications to the rapid test or to the test protocol (including sampling) may only be made following advance notification to the EURL, and provided that the EURL finds that the modification does not reduce the sensitivity, specificity or reliability of the rapid test. That finding shall be communicated to the EU Commission and to the national reference laboratories (based on Regulation (EC) No. 956/2010 amending Regulation (EC) No 999/2001).

## Description and Principles

This kit uses a proprietary method licensed to Microsens Biotechnologies (London, U.K.; patent pending) that allows detection of abnormal prions. A PrP<sup>Sc</sup>-specific ligand is immobilized on the surface of the antigen-capture plate. Test samples are prepared by homogenizing the tissues and then diluting the sample with working plate diluent. After the sample is applied to the plate, the disease-associated conformer binds to the immobilized ligand with high affinity. The plates are washed to remove unbound materials, including the normal conformer of the PrP protein. Following incubation with conditioning buffer, the captured antigen is then detected using a PrP-specific antibody that has been conjugated to horseradish peroxidase (HRPO). The plate is washed to remove unbound conjugate and a peroxidase substrate is added. Color development is related to the relative amounts of PrP<sup>Sc</sup> captured by the ligand immobilized in the microtiter plate well.

Interpretation of sample results is based on the sample absorbance. A sample whose  $A_{450} - A_{REF}$  is less than the cutoff value is considered to be negative by the IDEXX HerdChek BSE-Scrapie Antigen Test Kit. Samples whose  $A_{450} - A_{REF}$  is greater than or equal to the cutoff are classified as positive for PrP<sup>Sc</sup>. A confirmatory assay such as immunohistochemistry is required for all positive test results.

## Kit Components

Store all components at 2–8°C.

| Item  | 460 tests  |
|---|------------|
| <b>A</b> Antigen-capture plates   | 5 plates   |
| <b>N</b> Negative control—Nonreactive with antigen-capture plate; preserved with sodium azide             | 5 x 1 mL   |
| <b>P</b> Positive control—Noninfectious, reactive with antigen-capture plate; preserved with sodium azide | 5 x 1 mL   |
| <b>D1</b> Plate diluent component 1; preserved with sodium azide  | 20 mL      |
| <b>D2</b> Plate diluent component 2   | 5 x 200 µL |
| <b>R</b> Reconstitution diluent   | 20 mL      |
| <b>CB</b> Conditioning buffer; preserved with sodium azide  | 60 mL      |
| <b>CC</b> Conjugate concentrate; preserved with Bronidox L and methylisothiazolone                        | 300 µL     |

|               |   |             |
|---------------|---|-------------|
| <b>SRB-CC</b> | Small-ruminant brain conjugate concentrate; preserved with Bronidox L and methylisothiazolone | 300 $\mu$ L |
| <b>CD</b>     | Conjugate diluent buffer with detergents and protein stabilizers; preserved with Proclin      | 60 mL       |
| <b>W1</b>     | 10X wash solution 1; preserved with sodium azide  | 450 mL      |
| <b>W2</b>     | 10X wash solution 2; preserved with gentamicin  | 450 mL      |
| <b>T</b>      | TMB substrate   | 60 mL       |

**Notes:** See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

See the end of this insert for reagent hazard and precaution warnings.

### Materials and Equipment Required (Not Provided in Kit)

- Precision pipettes and multichannel pipettes suitable for delivering between 25 and 200  $\mu$ L. Reagent volumes listed in the Test Procedure require pipette precision less than or equal to 5%.
- Graduated cylinders for wash solutions
- Hard plastic or adhesive plate covers and reagent-dispensing trays
- Disposable sample collection dissection instruments or sample collection device
- 96-well plate reader (equipped with 450-nm and 620–650-nm filters) and washer
- FastPrep\* (FP120A, FP220A), FastPrep\*-24, Precess 48\*, Precess 24\* or Precellys 24\* instrument.
- Accessory pack—dilution plates, tissue-disruption tubes with beads, extended-length homogenate transfer pipette tips and adhesive plate covers (available from IDEXX)
- Personal protective equipment—safety glasses, lab coats, disposable gloves, shoe covers, hairnets and face masks
- 0.5–1.0 N HCl or 1.0 N H<sub>2</sub>SO<sub>4</sub> stop solution
- Sodium hypochlorite (bleach), 1N NaOH and 1N HCl, deionized water
- Optional—robotic sample processor with pipetting precision of  $\leq$ 2.5% as measured by CV analysis (Tecan, for example)
- Optional—microtiter plate shaker (e.g., IKA MTS 2/4)
- Optional—plate incubator that can maintain temperature of 32–37°C and has minimal air flow
- Optional—dry heat block unit for 1.5–2 mL microfuge tubes (capable of maintaining 70°C temperature)
- Optional—1.5–2 mL unskirted conical screw-cap microfuge tubes


### Tissue Sampling and Preparation

#### A. Bovine and Small Ruminant Brain Tissue

1. Sampling and laboratory testing must follow the regulation (EC) No 999/2001, Annex X, Chapter C, which refers in terms of collection of samples to the latest edition of the *Manual of Standards for Diagnostic Tests and Vaccines* of the International Office of Epizootic Diseases (OIE) stating: “The preferred sample for immunoassay should be at or as close to the obex as possible, but no further than 1.5 cm anterior to the obex.” Collect 0.30 g ( $\pm$ 0.05 g) of nervous tissue from the left or right side of the brain stem from the level of the obex (whenever possible) with dissection instruments, and weigh the sample to ensure the correct amount. Alternatively, the IDEXX sample collection device can be used for bovine obex collection as described in the Appendix. Personnel collecting the obex must be trained in the sampling method.

The diagram to the right indicates the correct sampling area.

**NOTE:** After sample collection, a complete hemisection of the brain stem with an intact obex region and cerebellum (if present) must remain available for confirmatory testing.

**Figure 1.** Bovine brain stem.

© British Crown Copyright (February 2005) reproduced with kind permission of the Veterinary Laboratories Agency.

2. Place the tissue in a tissue-disruption tube and cap tightly. Tubes are provided with ceramic beads and buffer.
3. Four tissue-disruption instruments have been validated for use with the IDEXX BSE-Scrapie EIA. Place the tubes into the instrument and disrupt as indicated for the applicable instrument. If grinding is insufficient, repeat one cycle.
  - FastPrep\* instrument program: Grind the samples for 40 seconds at the maximum speed (6.5 m/s). If a second cycle is required, the instrument must be allowed to cool 5–10 minutes between cycles.
  - Precess 48\*, Precess 24\* and Precellys 24\* instrument program: Grind the samples 2 times for 20 to 25 seconds at 6500rpm, with a 5 second delay between cycles.
4. Homogenates (fresh or thawed) can be held at 18–26°C for up to 4 hours prior to the start of the assay.

The number of samples to be prepared in a single session is flexible. Homogenates can be stored for up to 24 hours at 2–8°C or held at ≤-20°C for up to 6 months. Frozen homogenates must be thawed and thoroughly mixed by inversion before use. Tissue samples can be stored at -80°C.

### **B. Small Ruminant Lymph Node or Spleen Tissue**

1. Collect 0.30 g (±0.05 g) of tissue from lymph node or spleen tissue. For lymph nodes, tissue should be sampled to maximize the presence of germinal center cells. Mince the tissue into 8–10 small pieces.
2. Proceed with processing and sample storage as described for brain tissue.  
**NOTE:** Lymph node or Spleen Tissue cannot be used in the context of official sampling and testing in the framework of Regulation (EC) No 999/2001.

### **C. Cervid Brain and Lymph Node Tissue**

1. Prepare the tissue disruption tubes by removing 0.2 ml of solution from each tube.
2. Collect 0.30 g (±0.05 g) of tissue from retropharyngeal lymph node or brain (obex preferred). Brain and lymph node tissues should be tested individually. For lymph nodes, tissue should be sampled to maximize the presence of germinal center cells. Mince lymph node tissue into 8–10 small pieces.

Please note that in Germany individual testing of both brain and lymph node tissue must be undertaken. This decision follows the EURL guidelines (TSE EU Reference Laboratory Guidelines for the detection of Chronic Wasting disease in cervids, <https://science.vla.gov.uk/tse-lab-net/documents/cwd-guidelines.pdf> and Scientific opinion on chronic wasting disease (II), EFSA Panel on Biological Hazards (BIOHAZ), EFSA Journal 2018;16(1):5132. 59 pp.).

3. Place the tissue in a tissue-disruption tube and cap tightly.
4. Place the tubes into the tissue disruption instrument and disrupt as indicated for the applicable instrument. If grinding is insufficient, repeat one cycle.

FastPrep\* instrument program: Grind the samples for 30 seconds at the maximum speed (6.5 m/s). Rest the instrument for five minutes. Repeat grinding for another 30 seconds at maximum speed. If grinding is insufficient by visual inspection, repeat one cycle.

Precess 48\*, Precess 24\* and Precellys 24\* instrument program: Grind the samples for 30 seconds at 6500rpm. Rest the instrument for five minutes. Repeat grinding for another 30 seconds at 6500rpm. If grinding is insufficient by visual inspection, repeat one cycle.

## Preparation of Reagents

---

### Wash Solutions (Wash Solution 1, Wash Solution 2)

The wash solution concentrates should be brought to 18–26°C and mixed to ensure dissolution of any precipitated salts. Each wash concentrate must be diluted 1:10 with distilled or deionized water before use (e.g., 40 mL of concentrate plus 360 mL of water per plate to be assayed).

### Plate Diluent Component 2

Plate diluent component 2 (D2) is supplied as a lyophilized compound. The solution is prepared by adding 200  $\mu$ L of reconstitution diluent (R), letting it stand for approximately 1 minute and then mixing gently. Use within 1 hour of preparation.

### Working Plate Diluent

Plate diluent component 1 (D1) should be brought to 18–26°C. Prepare the working plate diluent by adding 1 part plate diluent component 2 (D2; prepared as above) to 25 parts plate diluent component 1 (D1) and mix thoroughly (e.g., 120  $\mu$ L D2 to 3.0 mL D1 per plate). Approximately 2.75 mL of working plate diluent is required per plate when testing brain samples. Lymph node or spleen testing requires approximately 5 mL of working plate diluent per plate. Working plate diluent should be prepared and used the same day.

### Negative and Positive Controls

Negative and positive controls are supplied lyophilized. Reconstitute each control by adding 1 mL of reconstitution diluent. Allow the solution to stand for approximately 1 minute and then mix vigorously. Use within 2 hours of preparation. **DO NOT DILUTE NEGATIVE OR POSITIVE CONTROLS INTO WORKING PLATE DILUENT.**

### HRPO-Conjugated Anti-PrP Antibody Solutions

The HRPO-conjugated anti-PrP antibody solutions are prepared by diluting the appropriate conjugate concentrate (see note below) into the conjugate diluent (CD) as indicated on the label (for example, a 1:100 dilution would require 120  $\mu$ L of conjugate concentrate to 12 mL of conjugate diluent). Refer to the conjugate concentrate label for the correct dilution factor. Diluted conjugate should be prepared and used within 4 hours.

**IMPORTANT NOTE:** There are two conjugate concentrates available for use with this test. Select the appropriate one for the tissue to be tested:

- **Conjugate concentrate (CC)**—Use this conjugate when testing bovine and cervid brain samples, small-ruminant and cervid lymph node samples and small ruminant spleen samples.
- **Small-ruminant brain conjugate concentrate (SRB-CC)**—Use this conjugate when testing Small ruminant brain tissue.
- Negative and positive control wells must be included for each type of conjugate tested.

### Acid Stop Solution

Assay stop solution is not provided in the test kit. Stop solution (0.5–1.0 N HCl or 1.0 N H<sub>2</sub>SO<sub>4</sub>) can be purchased at working concentration or prepared from concentrate.

All three protocols described below require that reagents be at 18–26°C before use. Before starting the test, prepare the solutions to be used in the assay. Mix all reagents by gently swirling. Controls (negative and positive) should be mixed vigorously and tested in duplicate. A plate cover should be used to cover the plate for the duration of the assay.

## Storage of Prepared Reagents

---

| Item                          | Reconstitution Volume | Shelf Life                                |
|-------------------------------|-----------------------|---|
| N/P Negative/Positive control | 1 mL                  | 2 hours at 18–26°C<br>(6 months at -20°C) |
| D2 Plate diluent 2            | 200 µL                | 1 hour at 18–26°C<br>(6 months at -20°C)  |
| Working plate diluent         | NA                    | 8 hours at 18–26°C                        |
| HRPO:anti-PrP solutions       | NA                    | 4 hours at 18–26°C                        |
| Wash solution 1–1X            | NA                    | 1 week at 18–26°C                         |
| Wash solution 2–1X            | NA                    | 1 week at 18–26°C                         |

**Store any unused portion of plates in a dark, desiccated, sealed container.**

## Test Procedure

---

Sample homogenates are prepared as described in the Tissue Sampling and Preparation section. **A robotic sample processor can be used in place of the manual method from Step 1 or once the controls and diluted samples have been added to the antigen-capture plate (Step 3).**

**Important: Cover each assay plate with a solid plastic or adhesive plate cover during all reagent incubations. If reagent incubations are conducted in a biosafety cabinet, plates must be covered with adhesive sheets.**

### Assay Protocols

The IDEXX BSE-Scrapie EIA has two approved protocols for brain tissue: Short and Ultra-Short. The protocols have equivalent performance but have varying equipment requirements for decreased assay time. The protocols are detailed in the table on the next page.

**NOTE:** The protocol for small-ruminant spleen and lymph node testing is different from the Short and Ultra-Short protocols and is described in the Assay Protocol table below.

### Dilution of Sample in Working Plate Diluent

Set up a template indicating where the sample positions are located on the antigen-capture plate and the dilution plate. Reserve duplicate wells for the kit controls. Working plate diluent can be added to the dilution plate before or after the sample. The ratio is 30 µL of working plate diluent per 120 µL of sample homogenate (small-ruminant spleen and lymph node ratio is 50 µL diluent per 100 µL sample homogenate).

Resuspend the homogenates by inversion and then carefully pipette the sample by extending a homogenate transfer pipette tip through the beads and drawing sample from the tissue-disruption tube. Carefully dispense each sample into the dilution plate, taking care to avoid creating bubbles in the homogenate or leaving any residual homogenate on the edges of the dilution plate well.

After the homogenate is diluted, thoroughly mix samples, taking care to avoid bubbles. Mixing can be done with a pipette or on a plate shaker. If a shaker is used, it will be necessary to optimize the speed and time to ensure complete mixing without splashing of the samples. Proceed with the assay within 2 hours.

### Use of Plate Shaker for Sample Incubation (Short and Ultra-Short Protocols)

The Short and Ultra-short protocols involve a slow rotation (200 ± 100 rpm) on a platform plate shaker for the sample incubation step only. The plate shaker should create a gentle horizontal, circular motion. Although the sample within each microtiter plate well will move slightly, it may not be visually apparent. The motion should not be vigorous enough to move sample up the side of the wells. Incubation times are decreased accordingly for sample and conjugate incubations as detailed in the following table. The Small Ruminant Spleen /Lymph node protocol requires that all incubations are stationary.

## Assay Protocols

| Test Procedure |  | Short Protocol   | Ultra-Short Protocol  | Small Ruminant Spleen/ Lymph Node Protocol                              |
|----------------|--|--|---|---|
| Sample Type    |  | Bovine, small ruminant and cervid brains, cervid lymph nodes   | Bovine and small ruminant brains  | Small ruminant lymph nodes and spleen                                   |
| Step           | Activity                                 | 18–26°C<br>All steps   | 32–37°C <sup>1</sup> : Incubations only (steps 3,5,8,10)<br>18–26°C: All other steps including wash steps | 18–26°C All steps   |
| 1              | Sample Addition to Dilution Plate        | 120 µL sample with 30 µL working plate diluent   |   | 100 µL sample with 50 µL working plate diluent                          |
|                |  | MIX WELL- see “Dilution of Sample in Working Plate Diluent” section.   |   |   |
| 2              | Sample Addition to Antigen Capture Plate | Pipette 100 µL of diluted sample onto the assay plate; mix controls, add 100 µL in duplicate; cover the plate with a plate cover.  |   |   |
| 3              | Capture Plate Incubation                 | 45–60 min.; slow shake 200 ±100 rpm  | 20–25 min.; slow shake 200 ±100 rpm   | 2–3 hours stationary  |
| 4              | Wash with 1X Wash 1                      | Wash the wells 6 times with ~350 µL 1X Wash 1.   |   |   |
| 5              | Conditioner: Addition/ Incubation        | Add 100 µL conditioning buffer to each well; cover the plate; incubate 10 ±1 min.  |   |   |
| 6              | Wash with 1X Wash 2                      | Wash the wells 3 times with ~350 µL 1X Wash 2.   |   |   |
| 7              | Conjugate Addition                       | Add 100 µL diluted conjugate (bovine brain and cervid samples use CC; small ruminant brain use SRB-CC) and cover the plate with a plate cover.   |   | Add 100 µL diluted CC conjugate and cover the plate with a plate cover. |
| 8              | Conjugate Incubation                     | 45–50 min.   | 25–30 min.  | 60–75 min.  |
| 9              | Wash with 1X Wash 2                      | Wash the wells 5 times with ~350 µL 1X Wash 2.   |   |   |
| 10             | Substrate Addition/ Incubation           | Add 100 µL substrate to each well; cover the plate; incubate 15 ±1 min. away from direct sunlight (do not use adhesive cover).   |   |   |
| 11             | Stop Addition/Read Plate                 | Add 100 µL acid stop solution; the plate can be held up to 30 min. in the dark prior to reading the optical density (450 nm) using a reference wavelength of (A <sub>REF</sub> ) 620 nm to 650 nm. |   |   |

1. Incubation at 32–37°C means placing the assay plate into an incubator that is pre-heated to 32–37°C.

## Overview of Test Formats

|  | Bovine Brain                   | Small-Ruminant Brain           | Small-Ruminant Lymph node and Spleen          | Cervid Brain & RPLN      |
|--|--------------------------------|--------------------------------|---|--------------------------|
| <b>Sample Prep:</b>  |                                |                                |   |                          |
| <b>Syringe Tissue Sampling</b>                                       | Yes                            | No                             | No  | No                       |
| <b>Tissue Mincing</b><br>(prehomogenization)                         | No                             | No                             | Yes   | Yes (RPLN)<br>No (Brain) |
| <b>Assay Conditions:</b>   |                                |                                |   |                          |
| <b>Conjugate Concentrate</b>   | CC                             | SRB-CC                         | CC  | CC                       |
| <b>Working Plate Diluent</b><br>(working plate diluent/sample ratio) | 30 $\mu$ L/120 $\mu$ L         | 30 $\mu$ L/120 $\mu$ L         | 50 $\mu$ L/100 $\mu$ L                        | 30 $\mu$ L/120 $\mu$ L   |
| <b>Approved Protocols</b>  | <b>Short &amp; Ultra Short</b> | <b>Short &amp; Ultra-Short</b> | <b>Small Ruminant Spleen &amp; Lymph Node</b> | <b>Short</b>             |
| <b>Test Cutoff</b>   | NC $\bar{x}$ + 0.120           | NC $\bar{x}$ + 0.180           | NC $\bar{x}$ + 0.180                          | NC $\bar{x}$ + 0.200     |
| <b>Optional Heat Treatment Protocol applies</b>                      | Yes                            | No                             | No  | No                       |

## Interpretation of Results

For the assay to be valid, the negative control mean (NC $\bar{x}$ ) must run with an  $A_{450}-A_{REF}$  value less than 0.150, and the positive control mean (PC $\bar{x}$ ) must have an  $A_{450}-A_{REF} \geq 0.400$ .

### Calculations

Calculation of negative control mean (NC $\bar{x}$ ):

$$NC\bar{x} = \frac{A1 (A_{450} - A_{REF}) + B1 (A_{450} - A_{REF})}{2}$$

Calculation of positive control mean (PC $\bar{x}$ ):

$$PC\bar{x} = \frac{C1 (A_{450} - A_{REF}) + D1 (A_{450} - A_{REF})}{2}$$

Calculation of the cutoff value:

Bovine cutoff = NC $\bar{x}$  + 0.120

Small-ruminant cutoff = NC $\bar{x}$  + 0.180

Cervid cutoff = NC $\bar{x}$  + 0.200

**NOTE:** See Step 11 of the Test Procedure Chart for the definition of  $A_{REF}$ .

### Results

Interpretation of sample results is based on absorbance for the sample. Samples with  $A_{450}-A_{REF}$  values less than the cutoff value are considered negative by the IDEXX HerdChek BSE-Scrapie Antigen Test Kit. Samples with  $A_{450}-A_{REF}$  values greater than or equal to the cutoff are classified as initially reactive for PrP<sup>Sc</sup>, and the homogenate should be retested in duplicate on the IDEXX HerdChek BSE-Scrapie Antigen Test Kit.

Bovine brain retesting may be done from the original tissue homogenate or from homogenate prepared using the optional bovine heat treatment protocol, described below. If either retest value is equal to or greater than the test cutoff, the sample is considered positive. The sample is considered negative when both retest replicates are less than the test cutoff. Within the EU all initially reactive bovine samples which are negative following the retest heat treatment protocol should be reported as negative.

Small-ruminant and cervid samples should be retested in duplicate directly from the original tissue homogenate - Do NOT heat treat small-ruminant or cervid homogenates. If either retest value is equal to or greater than the test cutoff, the sample is considered positive. The sample is considered negative when both retest replicates are less than the test cutoff.

In EU member states any rapid test positive samples must be sent to your countries National Reference Laboratory for TSEs or the European Union Reference Laboratory for confirmatory testing.

### **Optional Heat Treatment Protocol for Initially Reactive Bovine Homogenates: (The heat treatment protocol only applies to bovine samples.)**

Remove  $230 \pm 20 \mu\text{L}$  of the initially reactive tissue homogenate and dispense into a 1.5–2 mL conical bottom screw-cap vial. Place the vial in a dry heating block that has been preheated to  $70 \pm 2^\circ\text{C}$ . Heat the tube for  $10 \pm 1$  min. and then place the tube in an open rack at  $18\text{--}26^\circ\text{C}$  for a minimum of 20 minutes to allow the sample to cool. The homogenate should be retested within two hours of heat treatment, in duplicate, on the IDEXX HerdChek BSE-Scrapie Antigen Test Kit.

### **Precautions**

---

- Do not expose the tetramethylbenzidine (TMB) substrate to strong light or oxidizing agents. Use clean or disposable plasticware for dispensing TMB.
- Care should be taken to avoid contamination of kit components. Do not use components past their expiration dates, and do not intermix components from different kit lots.
- Some kit components contain sodium azide as a preservative (see the description of the kit components). Care should be taken to prevent contaminating the anti-PrP-HRPO conjugate with azide-containing solutions.
- Store all reagents at  $2\text{--}8^\circ\text{C}$ . Bring reagents to  $18\text{--}26^\circ\text{C}$  prior to use, and return to proper storage temperatures after use (see Storage of Prepared Reagents).
- Use separate dispensing trays for each reagent used in the assay. Avoid cross-contamination of the TMB substrate with the diluted conjugate solution. Do not pour unused TMB solution back into the bottle.
- Do not allow microtiter plates to sit more than 5 minutes between wash steps and the addition of reagents.

### **Safety Information**

---

- All personnel must receive appropriate training concerning risks related to BSE and scrapie, and recommended decontamination procedures. Biosafety procedures should be strictly followed as recommended by national safety regulations.
- The conditioning buffer contains chaotropic agents; avoid contact with skin and mucous membranes.
- The TMB substrate may irritate skin and eyes; avoid direct contact.
- Plate diluent 1 contains high concentrations of detergents; avoid direct contact.
- Avoid the use of glass containers in the lab.

### **Appendix**


---

#### **Bovine Obex Sampling with the IDEXX Sample Collection Device**

IDEXX has a sample collection device available as an alternative bovine obex extraction method. This device is a sampling syringe. The IDEXX sample collection device has been approved by the EURL. The guidance in this sampling section does not prohibit further information or instructions that are in line with regulation (EC) 999/2001 and amending regulations. When it is not possible to identify the correct anatomical area for sampling, use dissection instruments as described in the **Tissue Sampling and Preparation** section.

1. The brain stem should be collected at the abattoir via the occipital foramen using an appropriate tool or sample collection spoon. Identify the obex region of the brain as indicated by the presence of a V-shaped indentation on the upper surface of the brain stem (see the diagram in the Tissue Sampling and Preparation section). When it is not possible to identify the correct anatomical area for sampling, use dissection instruments as described in the Tissue Sampling and Preparation section.

2. Position the brain stem with the V-shaped section of the tissue facing upright. Place the tip of the sample collection syringe into the cut end of the caudal brain stem on the side that is to be sampled approximately 3 mm (enough so that it is securely lodged). It may be necessary to trim excess spinal cord tissue if the length from the base of the spinal cord to the apex of the V-shaped section is longer than 3–4 cm.

**Figure 2.** Cross section of the bovine brain stem at the level of the obex identifying the key target sites for tissue sampling: 1) solitary tract, 2) nucleus of the trigeminal nerve, 3) dorsal motor nucleus of the vagus nerve (diagram from *OIE Manual for Diagnostic Tests and Vaccines*, chapter 2.4.6.)

3. Hold the plunger of the syringe firmly. With your index finger, push the outer cylinder of the syringe into the brain stem, taking care not to allow the plunger portion to move in any direction. Refer to Figure 2 for proper alignment of the syringe to target the key tissue-sampling sites of the obex. As the syringe cylinder is pushed into the sample, it must stay within the selected side of the brain stem to prevent damage to the opposite side. A complete hemisection of brain stem with intact obex must be available for confirmatory testing.

4. The cylinder of the syringe will move through the brain stem and into the obex region. Ensure that the syringe has moved into the upper portion of the sampling area (see Figure 1). The cylinder should now contain the obex sample.

**NOTE:** The sample that you want (i.e., the obex area) is at the tip of the barrel.

5. Twist the cylinder to isolate the sample and carefully remove the syringe from the tissue.

6. If a significant portion of tissue is hanging from the end of the syringe, draw it into the cylinder by withdrawing the plunger slightly. The syringe can now be manipulated to remove air at the tip and close any gaps between sample portions.

**Note:** The inside of the syringe has several regularly spaced “detents” or grooves that can be located by feel as the plunger is moved through the syringe. The space between the detents provides an accurate measurement of sample volume.

7. When you have acquired the tissue in the syringe, push the plunger to align it with the nearest detent. The sample should be without gaps between the plunger and the tip of the syringe. Some excess sample material may extend beyond the tip of the syringe.

8. Wipe any minor residual tissue on the flat surface of the weigh boat. Do not press the plunger during this process because the sample will be ejected or compressed; both are undesirable.

9. Hold the tissue-disruption tube vertically in one hand and the syringe in the other, with the tip of the syringe just inside the mouth of the tissue-disruption tube. Dispense a measured amount of obex tissue into the disruption tube by moving the plunger through one detent and stopping on the second detent. The volume between detents is 150  $\mu$ L; a total of 300  $\mu$ L is dispensed into the tube (equivalent to 0.30 g  $\pm$  0.05 g of tissue).

10. Cap the tube and proceed with sample homogenization.

Personnel collecting obex samples with the IDEXX sample collection device should be well-trained on the use of the device to ensure collection of the correct area of the brain stem. Each technician should monitor sampling accuracy by conducting periodic checks of sample weight. A program of corrective action should be instituted when results fall outside of defined acceptance criteria. The IDEXX sample collection device is single-use and should be discarded after each sample to prevent cross-contamination.

### Limitation of Liability

To the maximum extent permitted by law, under no circumstances will IDEXX or our licensors be liable to you or any other person for loss of profit or use, special, incidental, consequential, indirect, exemplary, punitive or multiple damages, including without limitation for loss of goodwill, data or equipment or for business interruption, arising out of the manufacture, sale, supply or use of our products or services or failure or delay in delivering such products or services.

**For technical assistance:**

**IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4890**

**IDEXX Europe Tel: +800 727 43399**

\*HerdChek is a trademark or a registered trademark of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries. All other company and product names and logos are property of their respective holders.

Patent information: [idexx.com/patents](http://idexx.com/patents).

© 2023 IDEXX Laboratories, Inc. All rights reserved.

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SDS Update – Review online documents.

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**Positive control** – Harmful if swallowed. Harmful in contact with skin. Causes severe skin burns and eye damage. Causes serious eye damage. May cause respiratory irritation. Very toxic to aquatic life. Very toxic to aquatic life with long lasting effects. Wash thoroughly after handling. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection/hearing protection. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTRE/doctor. Collect spillage

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**H318 / H412 / P273 / P280 / P305+P351+P338 / P310 / P501**

**Negative Control** – Causes serious eye damage. Harmful to aquatic life with long lasting effects. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician. Dispose of contents/container in accordance with local/regional/national/international regulations.

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**H319 / P264 / P280 / P305+P351+P338**

**Plate Diluent 1** – Causes serious eye irritation. Wash thoroughly after handling. Wear eye protection/face protection. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

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**Plate Diluent 2** – Causes skin irritation. Causes serious eye irritation. May cause allergy or asthma symptoms or breathing difficulties if inhaled. May cause respiratory irritation. Avoid breathing dust. Avoid breathing vapours. Wash thoroughly after handling. Wear eye protection/face protection. Wear protective gloves. Wear respiratory protection. IF ON SKIN: Wash with plenty of water. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTRE/doctor if you feel unwell. If eye irritation persists: Get medical advice/attention. If experiencing respiratory symptoms: Call a POISON CENTRE/doctor. Take off contaminated clothing and wash it before reuse.

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**Conditioning buffer** – Harmful if swallowed. Harmful in contact with skin. Causes severe skin burns and eye damage. Causes serious eye damage. Harmful to aquatic life with long lasting effects. Do Not breath dust/vapors. Wash thoroughly after handling. Do not eat, drink or smoke when using this product. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician.

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H317 / P261 / P272 / P280 / P333+P313 / P362+P364 / P501

**Conjugate concentrate / SRB conjugate concentrate** – May cause an allergic skin reaction. Avoid breathing mist/vapors. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. Dispose of contents/container in accordance with local/regional/national/international regulations.

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H317 / H412 / P261 / P273 / P280 / P302+P352 / P333+P313 / P362+P364 / P501

**Conjugate Diluent** – May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Avoid breathing mist/vapors. Avoid release to the environment. Wear protective gloves. IF ON SKIN: Wash with plenty of water. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. Dispose of contents/container in accordance with local/regional/national/international regulations.

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© British Crown Copyright (février 2005) reproduit avec l'amiable autorisation de la Veterinary Laboratories Agency.
Copyright © Corona británica (Febrero de 2005). Reproducida con el consentimiento de la Agencia de Laboratorios Veterinarios de GB (VLA)
© British Crown Copyright (Februar 2005) Wiedergabe mit freundlicher Genehmigung des Veterinary Laboratories Agency.
© Copyright British Crown (Febbraio 2005) riprodotto grazie al cortese permesso rilasciato dal Veterinary Laboratories Agency.
© British Crown Copyright (Fevereiro 2005), reproduzido com a autorização e cortesia da Agência de Laboratórios Veterinários.
IDEXX B.V.